

Luna EMG

clinical
applications

EGZOTech



The content presented in this document, including that relating to our products, is addressed to specialists – medical professionals and healthcare professionals.

Who is it for?

- Neurological patients
- Orthopedic patients
- Oncology patients
- Pediatric patients
- Bedridden patients



EGZOTech



The content presented in this document, including that relating to our products, is addressed to specialists – medical professionals and healthcare professionals.

Patient Qualification for Therapy

Indications

- ✓ **Neurology:** Patients with stroke, spinal cord injuries, multiple sclerosis (MS), and cerebral palsy (CP).
- ✓ **Orthopedics:** Conditions after surgery and injuries
- ✓ **Bed-site treatment:** Intensive Care Units (ICU), post traumatic, palliative care, geriatric patients immobilized in bed
- ✓ **Pediatrics:** Muscular dystrophies, motor neuron diseases, myopathies, palsies.

Contraindications

- ✗ **General contraindications:** osseous and joint instability, non-consolidated fractures, ataxia, unstable vital functions
- ✗ **Significant limitation of range of motion:** significant spasticity, fixed contractures, anatomical features that prevent the positioning and stabilization of the patient
- ✗ **Cognitive impairment :** lack of contact/cooperation preventing active participation and generation of movement intention

Example Clinical Conditions

Neurology

- Stroke
- Spinal cord injury
- Peripheral nerve injuries
- Demyelinating diseases
- Paresis

Bed-site treatment

- ICU-Acquired Weakness
- Critical Illness Myopathy & Neuropathy
- Severe deconditioning & prolonged bed rest
- Acute stroke - early bedside mobilization

Orthopedics

- Reconstructions (ACL/PCL, rotator cuff)
- Post-operative conditions
- Post-traumatic conditions
- Following limb amputation

Other

- Pediatrics (Cerebral Palsy, SMA, Dystrophies)
- Geriatrics & Sarcopenia
- Oncology & Cancer-Related Fatigue
- Rheumatoid & Inflammatory Joint Diseases
- Respiratory Muscle Training

EGZOTech



The content presented in this document, including that relating to our products, is addressed to specialists – medical professionals and healthcare professionals.

↑ What Are Customers Looking For?

tune

Universal system for assessment and therapy

accessibility_new

Engaging patient training

arrow_upward

A bridge between passive and active therapy

spa

Relieving the therapist of repetitive, hard work



EGZOTech



The content presented in this document, including that relating to our products, is addressed to specialists – medical professionals and healthcare professionals.

Clinical Application 1: Foot Drop

Stage 1: Flaccid Paralysis

✓ Patient Qualification / Clinical Profile

Complete absence of active motor response in ankle dorsiflexors. Inability to generate trace movement. Flaccid paralysis of lower extremity muscles with no spasticity.

Lovett: 0 (No contraction) | Ashworth: 0 (Flaccid)

→ Primary Rehabilitation Goals

Preserve passive ankle Range of Motion (ROM) to prevent joint contractures caused by un-opposed antagonist muscles. Activate venous calf muscle pump and stimulate circulation.

Prescribed Luna EMG Exercises & Protocols

CPM & Progressive CPM

- **Exercise Mode:** Continuous Passive Motion (CPM) & Progressive CPM.
- **Duration:** 5–10 minutes per session.
- **Objective:** Mobilize ankle joint tissues and prevent contracture formation.

sEMG Muscle Monitoring

- **Target Muscle:** Tibialis Anterior peroneus muscles
- **Monitoring:** Real-time sEMG tracking during passive movements to detect early neuromuscular firing or reflex activity.

Venous Pump Activation

- **Mechanism:** Alternating ankle dorsiflexion and plantarflexion.
- **Role:** Serves as a warm-up protocol prior to secondary therapeutic interventions.

Clinical Application 1: Foot Drop

Stage 2: Flaccid Paresis

✓ Patient Qualification / Clinical Profile

Trace muscular contraction or minimal non-functional active movement present in dorsiflectors. Inability to overcome limb gravity or achieve functional joint range.

Lovett: 1–2 (Trace / Slight) | Ashworth: 0

→ Primary Rehabilitation Goals

Promote neuromuscular re-education and assist active movement initiation. Maintain passive/active ROM, stimulate calf muscle pump, and enhance motor recruitment.

Prescribed Luna EMG Exercises & Protocols

Warm-up & Preparation

- **CPM / Progressive CPM:** 5–10 minutes warm-up to prepare tissue and mobilize joint before active biofeedback work.

sEMG Biofeedback & Games

- **Tibialis Anterior:** More focus on short, rapid bursts (<2s) with high threshold (>80% MVC), but also longer contractions
- **Peroneus Longus:** Phase 1 short bursts; Phase 2 sustained holds (10–15s @ 50–60% MVC).

Reactive EMG & Active Control

- **Reactive EMG (T&R / T&H):** Dual-channel mode (TA + PL as agonists) for dorsiflexion with lateral stability.
- **Proprioception:** Dynamic Reversal if active motion occurs.

Clinical Application 1: Foot Drop

Stage 3: Active Flaccid Paresis

✓ Patient Qualification / Clinical Profile

Active ankle dorsiflexion movement restored against gravity and moderate resistance without hypertonia. Voluntary control returning, requiring progressive strength loading.

Lovett: 3-4 (Fair to Good) | Ashworth: 0

→ Primary Rehabilitation Goals

Expand active dorsiflexion Range of Motion (ROM). Build muscle force, power, and endurance. Refine foot stability, dynamic direction changes, and ankle proprioception.

Prescribed Luna EMG Exercises & Protocols

Warm-up

- **CPM / Progressive CPM:** 5–10 min if joint stiffness is present prior to active loading.

Resisted EMG Biofeedback

- **Biofeedback Games:** Target threshold recruitment under added mechanical resistance.
- **Resistance setup:** Elastic band around feet during dorsiflexion and pronation.

Active & Resisted Strength

- **Dynamic Reversal:** Bi-directional ankle movement training.
- **Strength Training:** Elastic resistance & Weight Lifting tailored to patient capabilities.

Clinical Application 1: Foot Drop

Stage 4: Spastic Paresis Management

✓ Patient Qualification / Clinical Profile

Dorsiflexor muscle weakness combined with elevated hypertonia/spasticity in Gastrocnemius/Soleus complex. High resistance during passive movement requires strict precautions.

Lovett: 1-2 (Trace / Weak) | Ashworth: 1-2 (Low / Moderate)

→ Primary Rehabilitation Goals

Maintain/improve joint ROM while actively suppressing calf spasticity. Promote reciprocal inhibition between Tibialis Anterior and Gastrocnemius. Increase strength safely.

Prescribed Luna EMG Exercises & Protocols

Low-Speed CPM & Safety Guard

- **Low Speed CPM:** Smooth passive movement to avoid triggering spastic stretch reflexes.
- **Safety Threshold:** Enable "Return on Max Torque" (3s limit) to handle spastic surges.

Agonist / Antagonist Biofeedback

- **TA & Peroneus Training:** Active dorsiflexion and inversion recruitment.
- **Triceps Surae Training:** Voluntary muscle relaxation focus (15–40s hold @ 20–30% MVC).

Reactive EMG (Agonist & Antagonist)


- **Agonist Mode:** Dorsiflexion with lateral stability (TA + Peroneus Longus).
- **Antagonist Mode:** Alternating Tibialis Anterior activation with Gastrocnemius relaxation.

Clinical Application 1: Foot Drop

Electrodes Positioning:

- red - tibialis anterior,
- blue - peroneus longus
- yellow - peroneus brevis

reference - knee cap



EGZOTech



The content presented in this document, including that relating to our products, is addressed to specialists – medical professionals and healthcare professionals.

Clinical Application 2: Stroke

Stage 1: Early Stage (24–48h to 2 Weeks Post-Stroke)

Patient Qualification

Patients in the early post-stroke phase, immediately following vital signs stabilization.

Flaccid paresis

Vital stability

Initiation: Start as early as medically feasible.

Therapeutic Objectives

- Activate muscle pump and improve peripheral circulation.
- Prevent edema in distal limb extremities.
- Early neurorehabilitation and prevention of immobility complications.

Exercise Program

- **Bedside Exercises:**
 - Continuous Passive Motion (CPM)
 - Reactive EMG (T&R / T&H)
- **sEMG Isometrics:**
 - EMG Biofeedback, Advanced Biofeedback, EMG Games
 - Short isometric contractions (tailored to prevent fatigue)
 - Optional sEMG Preview without enforced threshold

Clinical Application 2: Stroke

Stage 2: Subacute Stage (2 Weeks to 6 Months)

Patient Qualification

Window of high neuroplastic potential. Patients exhibiting gradual return of motor control and strength.

Spasticity Risk

Associated Movements

Precaution: Do not introduce resistance training too early to avoid triggering synkinesis!

Therapeutic Objectives

- Isolate movements and inhibit unwanted synkinesis.
- Promote correct movement patterns without compensation.
- Enhance muscle strength and neuromuscular control.

Exercise Program

- **2-Channel Reactive EMG (Isolation):**
 - *Relaxation / Antagonists:* Ankle dorsiflexion (Tibialis Anterior) + antagonist relaxation (Biceps Femoris / Gastrocnemius)
 - *Relaxation:* Arm flexion (Deltoid) + Trapezius relaxation
 - *Agonists:* Dorsiflexion with lateral stabilization (ankle)
- **sEMG Endurance:** Contractions up to 15s (50% MVC threshold).
- **Later Phase:** Elastic Resistance, Dynamic Reversal.

Dosage: 45 min (2-3 min breaks every 15 min), 5 days/week.

EGZOTech



The content presented in this document, including that relating to our products, is addressed to specialists – medical professionals and healthcare professionals.

Clinical Application 2: Stroke

Stage 3: Chronic Stage (> 6 Months Post-Stroke)

Patient Qualification

Chronic stroke patients (>6 months post-onset) with established movement habits and contractures.

Persistent Neuroplasticity

Fixed Compensations

Safety Rule: Work strictly within existing contracture ranges - never force joints!

Therapeutic Objectives

- Rebuild muscle mass and increase strength.
- Further eliminate compensatory motor patterns.
- Maximize functional independence in daily living.

Exercise Program

- **sEMG Biofeedback & Games:** Complex tasks demanding fine-grained muscle tension control.
- **Active & Resisted Training:**
 - Dynamic Reversal
 - Proprioception
 - Elastic Resistance and Weight Lifting
- **Ongoing Focus:** Modes suppressing synkinesis.

Model: Comprehensive (up to several hours daily). Robot as a supplement to training

Clinical Application 3: Total Knee Arthroplasty

Stage 1: Early Post-Op Stage (24h to 2 Weeks)

Patient Qualification

Early post-operative phase following total knee joint replacement surgery.

Acute Healing Phase

Swelling & Pain Focus

Early Mobilization: Essential for preventing adhesions and surgical complications.

Therapeutic Objectives

- Promote early mobilization and prevent post-operative complications (DVT, joint stiffness).
- Gradually improve passive and active joint Range of Motion (ROM).
- Activate distal muscle pump to improve peripheral venous circulation.
- Re-establish early motor unit recruitment in quadriceps and gluteal muscles.

Exercise Program

- **Knee CPM (Passive Motion):**
Tailored strictly to patient pain tolerance to preserve ROM.
- **Ankle CPM / Active Motion:**
Pumping exercises to activate distal leg muscles and boost circulation.
- **Advanced EMG Biofeedback Isometrics:**
 - Gluteal Muscle Isometrics: Short (<3s), medium (5–10s), and long (>15s) contractions.
 - Thigh Muscle Isometrics:
Performed in full knee extension with short, medium, and long contractions.

Dosage: Gentle sessions tailored to individual pain and swelling limits.

Clinical Application 3: Total Knee Arthroplasty

Stage 2: Subacute & Chronic Stage

Patient Qualification

Patients in the subacute and chronic phase post-TKA (2–8 weeks up to 6 months).

Healed Surgical Incision

Progressive Resistance Ready

Evidence-Based: Resisted strength training yields significantly superior recovery outcomes compared to unresisted active exercises.

Therapeutic Objectives

- Rebuild quadriceps and hamstring muscle mass and strength.
- Activate Vastus Medialis (VMO) in terminal extension range.
- Restore full functional knee extension and extension lag control.
- Enhance proprioception and neuromuscular control of the thigh.

Exercise Program

- **Knee CPM Warm-up:** 5–10 minutes passive motion before active training.
- **Reactive EMG (T&R / T&H):** VMO muscle activation specifically in the last 15° of knee extension.
- **Resisted Active Training:**
 - Dynamic Reversal, Elastic Resistance, Weight Lifting
 - Progressive load overload: 3–4 sets of 8–12 repetitions
- **Proprioception & Orthopedic Games:** Thigh muscle control exercises.

Frequency: Minimum 3 times per week with progressive overload.

EGZOTech



The content presented in this document, including that relating to our products, is addressed to specialists – medical professionals and healthcare professionals.

Quadriceps muscle - Electrode Placement



Bibliography

1. Thatte H, De Jesus O. Electrodiagnostic Evaluation of Peroneal Neuropathy. StatPearls Publishing, Jan 2023,
2. Vive S, Bunketorp-Käll L, Carlsson G. Experience of enriched rehabilitation in the chronic phase of stroke. *Disabil Rehabil.* 2022;44(3):412-419. doi: 10.1080/09638288.2020.1768598.
3. Doumen S, Sorba L, Feys P, Tedesco Triccas L. Efficacy and dose of rehabilitation approaches for severe upper limb impairments and disability during early acute and subacute stroke: a systematic review. *Phys Ther.* 2023;103(4):pzad002. doi: 10.1093/ptj/pzad002.
4. Liu Y, Yin JH, Lee JT, Peng GS, Yang FC. Early rehabilitation after acute stroke: the golden recovery period. *Acta Neurol Taiwan.* 2023;32(1):1-8.
5. Monsegue AP, Emans P, van Loon LJC, Verdijk LB. Resistance exercise training to improve post-operative rehabilitation in knee arthroplasty patients: A narrative review. *Eur J Sport Sci.* 2024;24(7):938-949. doi: 10.1002/ejsc.12114.
6. Konnyu KJ, Thoma LM, Cao W, Aaron RK, Panagiotou OA, Bhuma MR, et al. Rehabilitation for Total Knee Arthroplasty: A Systematic Review. *Am J Phys Med Rehabil.* 2023;102(1):19-33. doi: 10.1097/PHM.0000000000002008.
7. Tian J, Song M, Cho KJ, Lee HY, Ju SH, Lim JR, et al. Differences in Type 2 Fiber Composition in the Vastus Lateralis and Gluteus Maximus of Patients with Hip Fractures. *Endocrinol Metab (Seoul).* 2024;39(3):521-530. doi: 10.3803/EnM.2024.1935.
8. Lisi C, Caspani P, Bruggi M, Carlisi E, Scolè D, Benazzo F, Dalla Toffola E. Early rehabilitation after elective total knee arthroplasty. *Acta Biomed.* 2017;88(4-S):56-61. doi: 10.23750/abm.v88i4-S.5154.
9. Fransen M, Nairn L, Bridgett L, Crosbie J, March L, Parker D, Crawford R, Harmer AR. Post-Acute Rehabilitation After Total Knee Replacement: A Multicenter Randomized Clinical Trial Comparing Long-Term Outcomes. *Arthritis Care Res (Hoboken).* 2017;69(2):192–200. doi: 10.1002/acr.23117.